

01. BOOSTER PRO Guide

Taking the Lead Again with the All New 6-in-1 BOOSTER PRO

Revolutionary performer offering visible results in radiance, volume, elasticity, pores, and sonic vibration. All these aspects of care can be achieved with a single BOOSTER PRO.



Main Mode Overview

*Before use, be sure to familiarize yourself with the usage methods and recommended frequency for each mode.

Radiance Booster Mode

BOOSTER MODE

Usage Instructions Apply an appropriate amount of skincare product and attach the head to the skin. Roll gently in a circular motion to allow the product to be absorbed.

Recommended Frequency 1-3 times daily (5-10 minutes per session)

Line MC Mode

MC MODE

Usage Instructions Apply an appropriate amount of skincare product and attach the head to the skin. Use by moving from the inside of the face to the outside while keeping the head attached.

Recommended Frequency 1-3 times daily (5-10 minutes per session)

Elasticity Derma Shot Mode

DERMA SHOT MODE

Usage Instructions Apply an appropriate amount of skincare product and move slowly across the skin. Pause for 2 seconds at muscle points to relax the muscles.

Recommended Frequency 1-3 times daily (5-10 minutes per session)

Pores Air Shot Mode

AIR SHOT MODE

Usage Instructions After cleansing and drying the skin, gently glide the head's edge over the face.

Recommended Frequency 2-3 times a week (5 minutes per session)

Mode Sequence

*If different modes are being used on the same day, it is recommended to follow this sequence. Be sure to check the recommended frequency before use.

AIR SHOT MODE



BOOSTER MODE



MC MODE



DERMA SHOT MODE

Recommended Usage:
2-3 times a week

Recommended Usage: 1-3 times daily

LED INDICATION

*LED colors for each mode are preset, and LED color changes can only be made through the AGE-R app.



Blue



Orange



Green



Red




Purple

Mode-Specific Sonic Vibration

*The Sonic Vibration can only be turned ON/OFF through the AGE-R app.


Radiance Booster Mode

Vibrations are designed to help the skincare product absorb deeply.




Line MC Mode

Vibrations are synchronized with the lifting motion.




Elasticity Derma Shot Mode

Vibrations are suitable for muscle care.



Pores Air Shot Mode

Vibrations are intended for comprehensive pore care.



Skin Contact Sensor

*Due to the sensitive sensor, it may activate when there is a significant amount of skincare product on the head. For accurate skin recognition, please grasp the main part of the device as if surrounding it.



Smart and Safe Skin Contact Sensor Feature

Equipped with an advanced skin contact sensor that illuminates the LED light only when in contact with the skin, recording the duration of contact as the treatment time.

02. Bluetooth Connection

Achieve Full Control with AGE-R


From setting your preferred mode, intensity, volume, vibration, and LED colors to recording every aspect of your skincare journey with automatic reports. Use it effortlessly and keep track of your transformation.

Bluetooth Function Overview

01

Remote Control


Care made easy and convenient.



02

Automatic Reports


Turn your skincare routine into a habit.



03

Recording Before & After Use

Check for changes after 8 weeks of use.



* Usage records will be uploaded to the report only if the total usage time (skin contact time) for each mode is 20 seconds or more after turning the power on.

* The default volume setting is at level 3, and setting it to 0 (mute) will mute all sounds effects, including power on/off sound.

* The volume level and vibration ON/OFF status you change in the Control tab during use will be maintained as the default settings.



Get the most out of your BOOSTER PRO


Download the AGE-R app for free right now and experience the AGE-R BOOSTER PRO in a more professional and convenient way.

*Scan the QR code or search for 'AGE-R' in the app store.

Bluetooth Connection Instructions

*Make sure your smartphone's Bluetooth is turned on and the BOOSTER PRO is powered on to establish a connection.


01




02



03



04



- ① Access the AGE-R app and click on "Connect BOOSTER PRO."
- ② Turn on Bluetooth on your smartphone and power on the BOOSTER PRO.
- ③ After confirming the detected device, click on the "Connect" icon.
- ④ Verify that the "Connection Successful" screen and the Bluetooth connection sound from BOOSTER PRO are displayed correctly.

Precautions When Connecting via Bluetooth

1. Remote Control

- The Remote Control feature operates only when connected with the device, and functions such as LED color change, volume control, and vibration ON/OFF are available only within the Remote Control feature after Bluetooth connection.
- The volume level and vibration ON/OFF settings you configure in Remote Control will be saved as default values, but LED color changes will not be saved.

2. Automatic Reports

- Usage history when the device is used without being connected to the app is temporarily stored in the device and will automatically upload to the app's reports when connected.
- The intensity level in the report and mode-specific LED colors are saved as the longest used level and color.

3. Default LED Color Changes for Each Mode

- Changing the default LED color for each mode is possible only when the device is connected and can be done in the "Settings" tab.

4. Device Connection Failure

- Ensure that the device is not connected to another smartphone. If it is connected to another smartphone, disconnect it from the connected smartphone's app by going to [BOOSTER PRO Banner] → [Settings] → [Device Disconnect], and then try again.

5. Switching to a Different Smartphone

- If you want to connect the device to a different smartphone, first run the app on the currently connected smartphone. Go to [BOOSTER PRO Banner] → [Settings] → [Device Disconnect] to disconnect it, and then try connecting it to the new smartphone.



Only Available When Bluetooth is Paired

- Mode-Specific Color Changes
- Sonic Vibration ON/OFF
- Volume Control
- Checking Care Reports

03. Booster Mode

Booster Mode that Helps Enhance Radiance




BOOSTER MODE

Default Setting | Orange LED

Can be used with various types of skincare products and face masks. Please choose the most suitable method among the three touching methods below and use accordingly.

How to Use



Watering

Place the head against the skin and roll it in a circular motion, providing a boosting effect.



Boosting

Press and slide the head against the skin as if gliding, providing a boosting effect.



Focus

Place the head against the skin and press it on one spot for 1-2 seconds before moving to a next spot.

Precautions Before Use

*When using the product for the first time, start with the lowest setting and adjust as needed during use.



Prohibited Areas of Use

Do not use on the uvular area, open wounds, or areas where procedures or surgeries have been performed.



Recommended Areas for Low-Level Use

For sensitive areas such as the eyes, nose, and mouth area, it is recommended to use at a low level.



Booster Mode Usage Frequency

We recommend using it 1-3 times a day, following the recommended usage method.

Quick Start Guide



STEP 01

After cleansing, apply an appropriate amount of skincare product to the skin.



STEP 02

Ensure the head is in complete contact with the skin surface, then gently massage and absorb the skincare product.



STEP 03

Refer to the boosting method and continue using until the skincare product is fully absorbed. (1-3 times a day / 5-10 minutes per use)



STEP 04

Once the skincare product has been mostly absorbed, finish with targeted care on desired areas.

04. MC Mode

MC Mode for Line Volumizing Care



MC MODE

Default Setting | Green LED

Please use after applying the skincare product.

Adjust the amount of skincare product if you feel the stimulation to be too weak.

How to Use



LINING

Attach the head to the skin and use it in an upward motion.

TIP Enhancing the Effects of MC Mode

Please start with the MC Mode to help volumize fine and deeper lines. Then, switch to the Derma Shot Mode to massage the muscle for more elastic and rejuvenated skin.

MC MODE



DERMA SHOT MODE

It is recommended to use in this order.

Precautions Before Use

* When using the product for the first time, start with the lowest setting and adjust as needed during use.



Prohibited Areas of Use

Do not use on the uvular area, open wounds, or areas where procedures or surgeries have been performed.



Recommended Areas for Low-Level Use

For sensitive areas such as the eyes, nose, and mouth area, it is recommended to use at a low level.



MC Mode Usage Frequency

We recommend using it 1-3 times a day, following the recommended usage method.



STEP 01

After cleansing, apply an appropriate amount of skincare product to the skin.



STEP 02

Ensure the head is in complete contact with the skin surface, then move the head from the inside of the face to the outside in an upward motion.



STEP 03

For the neck and forehead areas, use it once more from bottom to top.



STEP 04

Once the skincare product has been mostly absorbed, finish with targeted care on desired areas.

05. Derma Shot Mode

Derma Shot Mode for Elasticity Care




DERMA SHOT MODE

Default Setting | Red LED

After applying skincare product, massage the area where muscle stimulation works best. While using it, you can observe the movement of your muscles directly.

How to Use



Focus

Slowly move the head in complete contact with the skin. Pause for 1-2 seconds at the points where muscle stimulation works best to relax the muscles.

TIP Enhancing the Effects of Derma Shot Mode

Please start with the MC Mode to help volumize fine lines and deeper lines. Then, switch to the Derma Shot Mode to stimulate the muscles for more elastic and rejuvenated skin.

MC MODE > **DERMA SHOT MODE** It is recommended to use in this order.

Precautions Before Use

* When using the product for the first time, start with the lowest setting and adjust as needed during use.



Prohibited Areas of Use

Do not use on the uvular area, open wounds, or areas where procedures or surgeries have been performed.



Recommended Areas for Low-Level Use

For sensitive areas such as the eyes, nose, and mouth area, it is recommended to use at a low level.



Derma Shot Mode Usage Frequency

We recommend using it 1-3 times a day, following the recommended usage method.



STEP 01

After cleansing, apply an appropriate amount of skincare product to the skin before use.



STEP 02

With the head in full contact with the skin, move slowly and pause for 1-2 seconds at areas where you feel the most stimulation, relaxing the muscles.



STEP 03

Massage gently from the neck line to the trapezius area and collarbone line.



STEP 04

Use for 5-10 minutes a day to care for a more elastic face.

06. Air Shot Mode

Air Shot Mode for Pore Care

AIR SHOT MODE

Default Setting | Blue LED

After cleansing, be sure to use it on dry skin, lightly sweeping it with the edge of the head. It won't properly transmit stimulation when used on wet skin.



How to Use



SWEEPING

Gently make contact with the skin using the edge of the head and sweep it in one direction.



For Broad Areas

When using on broad areas, sweep with the edge.



For Smaller Areas

For localized areas such as nose, chin, and between the eyebrows, gently sweep it with the tip.

Precautions Before Use

* When using the product for the first time, start with the lowest setting and adjust as needed during use.



Prohibited Areas of Use

Do not use on the uvular area, open wounds, or areas where procedures or surgeries have been performed.



Recommended Areas for Low-Level Use

For sensitive areas such as the eyes, nose, and mouth area, it is recommended to use at a low level.



Air Shot Mode Usage Frequency

We recommend using it once a day, 2-3 times a week, and taking a break for 2-3 days after each use.

Quick Start Guide



STEP 01

Use it on completely dry skin after cleansing.

**It won't function properly if the skin is wet.*



STEP 02

Gently make contact with the skin using the edge of the head and sweep it in one direction.



STEP 03

For localized areas such as nose, chin, and between the eyebrows, gently sweep it with the tip.



STEP 04

After 5 minutes of use, use the Booster Mode to absorb the skincare product for radiant skin.

07. Product Detailed Information

LCD Screen Overview



Default Screen

① Bluetooth Indicator



Bluetooth Disconnected



Bluetooth Connected
(*Integrated with dedicated AGE-R app)

② Battery Indicator



Remaining Battery




Low Battery
(*Requires Charging)

③ Current Usage Mode

④ Intensity Level

⑤ Care Time



Screen While Awaiting Mode Change After 5-Minute Care



① Automatic Shutdown Countdown



Will automatically shut down after 10 seconds if no button operation is performed.

② Additional Care Guidance

MODE button and app operation can be used to change to a different mode.



Charging/Low Battery Screen

① Charging Screen

Displayed in green, allowing you to check the charging status.

(*Does not operate during charging)

② Low Battery Screen

Displayed in red when the battery is low, indicating that charging is required.

③ Fully Charged Screen



The entire ring is displayed in green.